

Refractive Index Detector

NEW

Shodex RI-201

Higher Performance

A new refractive index (RI) detector series is now available. RI-201 is a highly sensitive RI detector incorporating a three-chamber flow cell.

■ Features

- A novel optical system (three-chamber flow cell) provides twice or better sensitivity than that of our previous detectors.
- The two-step temperature control method significantly reduces drift due to the fluctuation of room temperature.
- The detection limit for saccharides is approximately 2ng.

Refractive Index Detector
Shodex RI-201H

NEW

reasonable price

■ Features

- Adopts the same optical system as that of RI-101. (2 chamber-type)
- Reasonable price.


Product Code	F4010105	F4010106
Model	RI-201	RI-201H
	Analysis	
Flow Cell Type	3 Chamber-type	2 Chamber-type
Measurement Method	Deflection Type	
Refractive Index Range	1.00~1.75	
Measurement Range	0.125~256 μ RIU	0.25~512 μ RIU
Drift *	0.1 μ RIU/h	0.2 μ RIU/h
Linearity Range	$\geq 300\mu$ RIU	$\geq 600\mu$ RIU
Noise **	≤ 1 nRIU	≤ 2.5 nRIU
Response	0.1, 0.25, 0.5, 1, 1.5, 2, 3, 6sec	
Auto Zero	Full auto zero	
Auto Zero Range	All range	
Off-set Range	5 μ RIU	10 μ RIU
Off-set Resolution	25nRIU	50nRIU
Integrator Output (Sensitivity)	DC 0~1V (4mV/ μ RIU, 16mV/ μ RIU)	DC 0~1V (2mV/ μ RIU, 8mV/ μ RIU)
Cell volume	8 μ L	
Flow Rate	(Usual)	0.2~3.0mL/min
	(Max.)	10mL/min (Solvent; Pure Water)
Maximum Back Pressure	50kPa	
Internal Volume	IN→Cell ; 80 μ L Cell→OUT ; 600 μ L All (Cell→OUT) ; 690 μ L	IN→Cell ; 60 μ L Cell→OUT ; 600 μ L All (Cell→OUT) ; 670 μ L
	Recorder Output 0~10mV/FS	
External Input	-	
External Output	(1) READY (Temperature Control) (2) LEAK (3) ERROR (ROM, RAM, PARAMETER, HOME-POSITION, OVER-HEAT, OPT.-BALANCE, INTENSITY)	
Temperature Control	OFF, 30~50°C (1°C Step), 77°C Temp. Fuse (Double Temperature Control)	
Communication Port	USB	
Operator Support Function	None	
Wetted Materials	Stainless Steel 316, Teflon, Quartz Glass	
Power Source, Power Consumption	AC100~240V \pm 10%, 50/60Hz, 150VA max	
Dimension, Weight	W260 x D400 x H200 (mm), ca. 12kg	
Accessories	Power Cable, Signal Code, Connector Tube, Fuse, Operation Manual	

*Pure water 1mL/min, PURGE OFF



**Pure water, response: 1.5sec

■ Principle of New Optical System

In our previous optical system, a measurement light passing through the flow cell is refracted only once. The new three-chamber flow cell allows the light to be refracted twice, thereby increasing sensitivity at least twice. This doubles the degree of polarization light and thus halves the noise and the influence of drift.



Application



Refractive Index Detector

Shodex RI-100 Series

The Shodex RI-100 series includes versatile and highly sensitive RI detectors equipped with a color liquid crystal display, an automatic start-up function, and a validation wizard, and can be used with various HPLC systems.



■ Features

- Equipped with a color liquid crystal display; the detector monitors the status of the chromatograph in real time.
- The automatic start-up function automates complicated operations such as the substituting eluent in the reference cell and checking the baseline stabilization.
- The validation wizard achieves component validation with ease.
- Optimization of the temperature control method shortens the time for stabilization after the start-up and improves baseline stability.
- The leak sensor automatically stops the pump in case of solvent leakage.
- External input and output terminals and RS232C communication ports permit advanced automation of the system.


Product Code	F4010101	F4010104	F4010102	
Model	RI-101	RI-104	RI-102	
	Analysis	Semi-micro	Preparative	
Flow Cell Type	2 Chamber-type			
Measuring Method	Deflection Type			
Refractive Index Range	1.00~1.75			
Measurement Range	0.25~512 μ RIU	0.25~512 μ RIU	2.5~5120 μ RIU	
Drift *	0.2 μ RIU/h	0.2 μ RIU/h	2 μ RIU/h	
Linearity Range	\geq 600 μ RIU	\geq 600 μ RIU	\geq 6000 μ RIU	
Noise **	\leq 2.5nRIU	\leq 5nRIU	\leq 25nRIU	
Response	0.1, 0.25, 0.5, 1, 1.5, 2, 3, 6sec			
Auto Zero	Full Auto Zero			
Auto Zero Range	All Range			
Off-set Range	0~500mV (Same as Integrator Output)			
Off-set Resolution	10mV (Same as Integrator Output)			
Integrator Output (Sensitivity)	DC 0~1V (2mV/ μ RIU, 8mV/ μ RIU)	DC 0~1V (2mV/ μ RIU, 8mV/ μ RIU)	DC 0~1V (0.2mV/ μ RIU, 0.8mV/ μ RIU)	
Cell Volume	8 μ L	2.5 μ L	8 μ L	
Flow Rate	(Usual)	0.2~3.0mL/min	0.2~1.0mL/min	1.0~50mL/min
	(Max.)	10mL/min (Solvent; Pure Water)	1.0mL/min (Solvent; Pure Water)	100mL/min (Solvent; Pure Water)
Maximum Back Pressure	50kPa			
Internal Volume	IN→Cell: ca. 60 μ L Cell→OUT: ca. 600 μ L All (Cell→OUT): ca. 670 μ L	IN→Cell: ca. 10 μ L Cell→OUT: ca. 355 μ L All (Cell→OUT): ca. 370 μ L	IN→Cell: ca. 120 μ L Cell→OUT: ca. 510 μ L All (Cell→OUT): ca. 640 μ L	
Recorder Output	0~10mV/FS			
External Input	Purge On/Off, Auto Zero, Marker			
External Output	(1) READY (Automatic Start-up) (2) LEAK (3) ERROR (OVER HEAT/LOW LIGHT INTENSITY/NULL GLASS HOME POSITION/LOST PARAMETERS/OPTICAL BALANCE) (Contact Capacity: DC24V 0.1A max.)			
Temperature Control	OFF, 30~50°C (1°C step), 77°C Temp. fuse			
Communication Port	RS232C			
Operational Support Functions	(1) Automatic Start-up (Start-up Sequence) (2) Span/Validation Guide (3) Real Time Baseline Monitor			
Wetted Materials	Stainless Steel 316, Teflon, Quarts glass			
Power Source, Power Consumption	AC100~240 \pm 10%, 50/60Hz, 150VA max			
Dimensions, Weight	W260 x D400 x H200 (mm), ca. 13kg			
Accessories	Power Cable, Signal Code, Connector Tube, Fuse, Operation Manual			

*Pure water 1mL/min, PURGE OFF

**Pure water, response: 1.5 sec

Eluent Degassing Device

DEGASSER® ERC-3000α Series ERC-3115α, 3215α, 3315α, 3415α



*[DEGASSER®] is the registered trademark of E.R.C. Co., Ltd.

■ Principle

The liquid to be degassed runs through a membrane made of special synthetic resin in a vacuum chamber. Dissolved gas in the liquid passes through the membrane and is eliminated from the liquid because of its small molecular size, higher mobility than that of the liquid, and higher affinity with the special synthetic resin membrane.

■ Features

- Pressure sensor and leak monitor assure high safety.
- Dead volume is minimized to 7mL/channel.
- Any pump can be used because differential pressure is low.
- Continuous degassing mode helps eluent preparation for high-sensitivity analysis.
- EMC and LVD compliance, with CE marking.

Product Code	Y4617007	Y4617000	Y4617002	Y4617004
Model	ERC-3115α	ERC-3215α	ERC-3315α	ERC-3415α
Solvent Channels	1 Channel	2 Channels	3 Channels	4 Channels
Degassing Capacity	When ion-exchanged water saturated with air at 25°C is put through at a flow rate of 3 mL/min, no bubble is observed at the outlet of the apparatus. (Measured dissolved oxygen level at the outlet: 2 ppm at flow rate of 3 mL/min)			
Internal Volume	7mL/Channel			
Max. Flow Rate	20mL/min for each Channel (Eluent: 25°C Pure water)			
External Output	An open connector signal is delivered to the external output signal terminal, when "PRES" and "LEAK" LEDs light.			
Dimensions	W71 x D310 x H136 (mm)			
Power Source	AC100V~AC240V, 50/60Hz			
Functions, Displays	<ul style="list-style-type: none"> • Power On/Off Display: "POWER" LED lights, when the power is supplied. • Status monitoring function: "READY" LED lights, when the internal pressure in the vacuum chamber is below a predetermined limit. • Pressure monitoring function: "PRES" LED lights, when the internal pressure in the vacuum chamber does not reach a predetermined level within a predetermined time. • Leak monitor: "LEAK" LED lights, when the liquid leaks in the apparatus. • Self cleaning: The vacuum line is cleaned by air suction. • Vacuum pump operation switching function: NORM.: Controlled operation in normal run CONT.: Continuous operation in case of high degree of degassing 			
Weight	ca. 5.0kg	ca. 5.1kg	ca. 5.4kg	ca. 5.5kg

* For use of special solvents, such as fluorinated organic solvents, please contact Shodex or our distributors near you.