


Toyopearl GigaCap® CM-650 Resins

(weak cation exchange chromatography resins)



Part Numbers:	21946 Toyopearl GigaCap CM-650M, 75µm, 100mL	21951 ToyoScreen GigaCap CM-650M, 1mL x 6 ea.
	21947 Toyopearl GigaCap CM-650M, 75µm, 250mL	21952 ToyoScreen GigaCap CM-650M, 5mL x 6 ea.
	21948 Toyopearl GigaCap CM-650M, 75µm, 1L	21400 ToyoScreen Column Holder
	21949 Toyopearl GigaCap CM-650M, 75µm, 5L	42194 ToyoScreen Column Holder with fittings
	21950 Toyopearl GigaCap CM-650M, 75µm, 50L	

Product Description

Toyopearl chromatographic resins are comprised of a methacrylic polymer incorporating high mechanical and chemical stability. Resins are available as non-functionalized "HW" series resins for size exclusion separations, and derivatized with surface chemistries for alternative modes of chromatography such as ion exchange, hydrophobic interaction or affinity separations.

Toyopearl GigaCap CM-650 chromatographic resins are designed for cation exchange chromatography. This chromatographic mode separates molecules on the basis of ionic interactions between the sample and the resin. The separation is usually accomplished in buffered aqueous solution with a gradient of increasing ionic strength. Alternatively, pH adjustment may be used for control of elution.

Formulation

72% (v/v) slurry in 20% (v/v) ethanol

This sheet contains the suggested operating conditions and specification ranges for Toyopearl GigaCap CM-650 resins. Lot-specific specifications are included in the Certificate of Analysis (C of A) shipped with the product.

A. OPERATING CONDITIONS

Pressure Drop Across the Column	max. 3 bar (recommended)
Operating Linear Flow Rate	normally 10 – 600 cm/hour (depending on particle size)
Cleaning in Place/Sanitization	0.5mol/L NaOH or 0.1mol/L HCl
Long Term Storage Conditions	20% (v/v) ethanol (shipping solvent)

B. SPECIFICATIONS

Appearance	White resin slurry which settles upon standing
Particle Size Distribution (>80% within range)	50 – 100µm for M-grade
Mean Pore Diameter (base resin)	1000Å
Ion Exchange Capacity	0.225 ± 0.055 meq/mL
Protein Adsorption Capacity	≥ 110 mg/mL (human IgG)
Bacterial Count	≤ 100 CFU/mL
Endotoxin Concentration	≤ 10.0 endotoxin units/mL
Resin Volume per Container	The indicated volume is the settled resin volume.
Foreign Substance	not observed, no visual evidence
Eluable Matter	≤ 0.2%
Shelf Life Stability	Estimated minimum 10 years

RDS1020_MAY2009