

TSK-GEL® G2000SW Products

Part Numbers:	05788, 7.5mm ID x 30cm, 10µm 05102, 7.5mm ID x 60cm, 10µm 06727, 21.5mm ID x 30cm, 10µm 05146, 21.5mm ID x 60cm, 13µm	05371, Guard Column for P/N 05788 & 05102 05758, Guard Column for P/N 06727 & 05146 06819, TSKgel SW top-off for 7.5mm ID columns, 1g wet gel
----------------------	--	---

This sheet contains the recommended operating conditions and the specifications for TSK-GEL G2000SW columns. Installation instructions and column care information are described in a separate Instruction Manual.

A. OPERATING CONDITIONS	
1. Shipping Solvent:	0.05% NaN ₃ and 0.1M Na ₂ SO ₄ in 0.1M phosphate buffer, pH 6.7
2. Max. Flow Rate:	1.2mL/min (7.5mm ID) 8.0mL/min (21.5mm ID) When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so as not to exceed the maximum pressure drop. When changing solvents, use a flow rate equal to 25% of the maximum flow rate.
3. Standard Flow Rate:	0.5 - 1.0mL/min (7.5mm ID) 3.0 - 6.0mL/min (21.5mm ID)
4. Max. Pressure:	1MPa (21.5mm ID X 30cm) 2MPa (7.5mm ID X 30cm, 21.5mm ID x 60cm) 4MPa (7.5mm ID X 60cm)
5. pH Range:	2.5 - 7.5
6. Salt Conc.:	< 0.5 Molar
7. Organic Conc.:	0 - 100% for aqueous soluble organic solvents. Make gradual solvent changes using a shallow gradient at low flow rate.
8. Temperature:	10 - 30°C, Reduce flow rate when operating below 10°C.
9. Cleaning Solvents:	(1) conc. salt solution at low pH, e.g. 0.5M Na ₂ SO ₄ , pH 2.7 (2) methanol or acetonitrile in low conc. aqueous buffer (3) buffered solution of urea or guanidine NOTE: Choose a cleaning solvent based on sample properties, e.g. use (1) to remove basic proteins, and (2) to remove hydrophobic proteins. Chaotropic agents can solvate strongly adsorbed proteins, e.g. via hydrogen bonding.
10. Storage:	Store the column in mobile phase containing 0.05% NaN ₃ or 20% ethanol when it will not be used the next day. For overnight storage flush the column with mobile phase at low flow rate. Prevent air from entering the column!
11. Column Protection:	The use of guard columns is recommended to prolong the life of the analytical column. Guard column life depends greatly on sample cleanliness. As a general rule, guard columns should be replaced after every 30-40 sample injections, when the peaks become excessively wide, or when the peaks show splitting.
12. TSKgel SW top-off gel:	Occasionally due to accident, sample, mobile phase, or operational variables, a depression can develop at the column or guard column inlet. Use TSKgel SW top-off gel for filling in such voids.
B. SPECIFICATIONS	
The performance of TSK-GEL G2000SW columns is tested under the conditions described in the Data Sheet. All columns have passed the following quality control specifications:	
1. Number of Theoretical Plates (N):	≥ 10,000 (30cm columns) ≥ 20,000 (60cm columns)
2. Asymmetry Factor (AF):	0.7 - 1.6

DS1005 Revised 18AUGUST2009

Tosoh Bioscience LLC
156 Keystone Drive
Montgomeryville, PA 18936
Phone: (215) 283-5000
Orders & Technical Service: (800) 366-4875
FAX: (215) 283-5035



MD Scientific is authorized distributor
in Denmark for Tosoh Bioscience
www.md-scientific.dk
info@md-scientific.dk
Tel. 7027 8565