

The chemagic QA Software

Improve your process safety!

chemagic Magnetic Separation Module I with chemagic QA Software

You own a chemagic MSM I but are still doing the documentation of the samples or the used chemagic Kit manually. However it is very easy to equip your chemagic MSM I with the latest version of the chemagic QA Software.



Improve the safety of the entire process by using the chemagic QA software, which tracks every single step.



chemagic Dispenser chemagic MSM I chemagic QA Software

The chemagic QA Software

- QA means Quality Assurance
 - Barcode reading of the entire process, including the used chemagic Kit, the samples/primary samples and every single isolation/protocol step
- Simple handling
 - Self explanatory software that guides you through the isolation process
- Data management
 - Connection to an existing software (e.g. LIMS) is possible at anytime



Your benefits

- **Real time tracking**
 - **Tracking the progress of the running protocol**
 - **Status protocol gives information about the remaining time to external step or the end of the isolation process**

- **Increased process safety through reduced failure risk**
 - **Human errors are excluded, sample mix up is impossible**

- **Speed**
 - **Preinstalled programs enable starting the chemagic MSM I more quickly**

- **Service Protocols**
 - **Integrated wash protocols facilitate the maintenance of the chemagic MSM I**

**Do you have any further questions
or
need an individual proposal?**

I am looking forward to your inquiry.



**MD Scientific ApS
Klokkerfaldet 23
8210 Aarhus V
Denmark
Phone: +45 70 27 85 65
email: service@md-scientific.dk
www.md-scientific.dk**