

Interchim Flash Phases

Reverse phase, Normal Phase, Ion Exchange, Mixed mode and Chiral Columns from 2g to 1.5kg

Product	Sorbent	Phase	Feature	End-Capping	Porosity Å	Surface Area m ² /g	Particle Size µm	Shape
puriFlash	Silica	C18	18% C	Yes	60	500	30 & 50	Spherical
puriFlash	Silica	C18U	16% C	No	60	500	30 & 50	Spherical
puriFlash	Silica	C8	11% C	Yes	60	500	30 & 50	Spherical
puriFlash	Silica	NH ₂	5% C	No	60	500	30 & 50	Spherical
puriFlash	Silica	OH	Diol	No	60	500	30 & 50	Spherical
puriFlash	Silica	Si	Bare silica	No	60	500	30 & 50	Spherical
Atoll	PSDVB	XC	High Capacity	No	NA	1500	70	Spherical
Atoll	PSDVB	XWP	High Capacity	No	Wide Pore	1200	90	Spherical
Atoll	NA	ATH	Hydrophilic	No	70	800	30 & 75	Spherical
Atoll	PSHEMA	AEV	Hydrophilic/Hydrophobic	No	70	800	30 & 75	Spherical
Atoll	PSDVB	ATL	Hydrophobic	No	70	800	100	Spherical
Recovery	Silica	REC	C18, High Capacity	NA	NA	NA	50	Spherical
Recovery	Silica	RSI	Si, High Capacity	No	NA	NA	50	Spherical
Upti-Clean	Silica	LCC	10% C	Yes	60	500	50	Spherical
Upti-Clean	Silica	RPAQ	C18, 100% Hydrophilic	No	NA	NA	75	Spherical
Upti-Clean	Silica	C18-S2F	High Flow	Yes	60	500	140	Spherical
Upti-Clean	Silica	C18U-S2F	High Flow	No	60	500	140	Spherical
Upti-Clean	Silica	C8U-S	9% C	No	60	500	50	Spherical
Upti-Clean	Silica	C8-S2F	High Flow	Yes	60	500	140	Spherical
Upti-Clean	Silica	C2	6% C	Yes	60	450	60	Granular
Upti-Clean	Silica	PH-S	Phenyl, 9% C	No	60	500	50	Spherical
Upti-Clean	Silica	CH	10.5% C	No	60	450	60	Granular
Upti-Clean	Silica	Si-S	Bare silica	No	60	500	50	Granular
Upti-Clean	Silica	CN-S	8% C	Yes	60	500	50	Spherical
Upti-Clean	Silica	SCX	Strong acid, 0.7 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	WCX	Weak acid, 0.22 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	DEAE	Weak base 0.33 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	SAX	Strong base 0.30 meq/g	No	60	450	60	Granular

Upti-Clean	Silica	MM1	RP/SCX, 0.09 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	MM2	RP/WCX, 0.1 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	MM3	RP/SAX, 0.14 meq/g	No	60	450	60	Granular
Upti-Clean	Silica	MM4	RP/WAX, 0.15 meq/g	No	60	450	60	Granular
Upti-Clean	Florisil	FL	Standard grade	No	NA	NA	200	Granular
Upti-Clean	Florisil	FLPR	Pesticide residue grade	No	NA	NA	200	Granular
Upti-Clean	Polyamide	P6	6	No	NA	NA	100	Granular
Upti-Clean	Amberlite	XAD2	Polystyrene	No	90	330	300-1200	Spherical
Upti-Clean	Amberlite	XAD4	Polystyrene	No	50	750	300-1200	Spherical
Upti-Clean	Amberlite	XAD7	Acrylic ester	No	80	450	300-1200	Spherical
Upti-Clean	Amberlite	XAD16	Polystyrene	No	100	825	300-1200	Spherical
Upti-Clean	Alumina	ALA	Acidic	No	NA	200	32-63	Granular
Upti-Clean	Alumina	ALN	Neutral	No	NA	200	32-63	Granular
Upti-Clean	Alumina	ALB	Basic	No	NA	200	32-63	Granular
Upti-Clean	Silica	BCD	B-Cyclodextrine	No	NA	NA	40	Granular
Upti-Clean	Carbon	CG	Graphitized	No	Non-Porous	100	NA	Granular
Upti-Clean	Carbon	CA	Activated	No	NA	NA	40	Granular